

TRACKERTM Plus

Kyphoplasty System


Tracker PlusTM Kyphoplasty system is an advanced total solution for vertebral compression fractures due to osteoporosis and cancer.

The system has specific design for reducing spinal fractures, create a void and correct vertebral body deformity.



Tracker Plus

Kyphoplasty System




Tracker Plus Kyphoplasty System consist of 10 gauge bone access tools with various type of tips and 400-psi inflatable balloon catheter in volume options to make good usability in the treatment of percutaneous vertebral agumentiation procedures, such as balloon kyphoplasty

Tracker Plus

Balloon Catheter

- Available in 10 gauge
- High Pressure up to 27 ATM / 400 PSI




Product Code	Balloon Length(mm)	Rated Burst (ATM/PSI)	Balloon Diameter	Volume
GSB-1110	10	27/400	14	3
GSB-1115	15	27/400	15	4
GSB-1120	20	27/400	16	5

Tracker Plus

Flex Balloon Catheter


- Simple, comprehensive curved kyphoplasty procedure.
- The balloon is inflated, deflated and removed.
- The curved needle stays in place for cement delivery.



Product Code	Balloon Length(mm)	Rated Burst (ATM/PSI)	Balloon Diameter	Volume
GSFB-1110	10	27/400	14	3
GSFB-1115	15	27/400	15	4
GSFB-1120	20	27/400	16	5

Tracker Plus Balloon Inflation Device


- 20ml Disposable Device with a range of 0 to 24.6 ATM (0 to 350 PSI).



Product Code	Operation Range	Accuracy
GSD-1000	0 ATM (0 PSI) to 24.6 ATM (350 PSI)	+/- 10% of full scale

Tracker Plus Cement Dispenser System

- One step percutaneous bone access tools.
- 10-gauge working cannula and multiple types of tips bevel, trocar and diamond and thread.



Product Code	Product Description	Total length (mm)	Remark
GCD-N1100	Bevel Needle	184	GCD-N1000 (Kit)
GCD-N1200	Trocar Needle	184	
GCD-N1300	Diamond Needle	184	
GCD-N1400	Thread Needle	184	Option
GCD-S1000	Spacer	245	Option
GCD-F1000	Cement Filler	235	

Tracker Plus

Flex Cement Dispenser

- Simple, comprehensive curved kyphoplasty procedure.
- A unipedicular approach is created and a curved needle is inserted into an ideal location.



Product Code	Product Description	Total length (mm)
GCD-FD1000	Flex Drill	235
GCD-FF1000	Flex Cement Filler	245

Tracker Plus

Bone Biopsy Kit

- Tracker Plus Biopsy Kit is for collecting biopsy samples



Product Code	Product Description	Kit components
GSBB-00	Bone Biopsy Kit	1 x Syringe 1 x Biopsy Tools 1 x Conical tube

Tracker Plus **Bone Cement Delivery System**

- Tracker Plus Cement Delivery System makes it easier to mix bone cement and deliver it more accurately.



Cement Mixing System

Product Code	Product Description	Kit components
GSMD-00	Cement Mixing System	5 x Syringe 1 x Cement Mixer Kit

- One-step working channel.
- Curved or straight balloon and cement injection options.
- Streamlined packaging for both ASCs and hospitals.



Cement Injection System

Product Code	Product Description	Kit components
GSIG-00	Cement Injection System	1 x Cement Injection Device 1 x Cement Mixer Kit

TRACKER™ Plus
Kyphoplasty System