

PRECISION SPINE

DAKOTA ACDF™

STANDALONE CERVICAL SYSTEM


REDEFINING RESULTS

The Dakota ACDF™ System simplifies standalone cervical procedures to enhance the fusion process. An integrated titanium plate and PEEK cage interbody and instrumentation facilitate intraoperative challenges and decrease OR time while providing reproducible results.

- **Dual thread screws** maximize cortical and cancellous bone interface
- **One-step, integrated locking mechanism** provides tactile and visual confirmations for increased confidence
- **Intuitive, low profile instrumentation** accommodates diverse patient anatomies



45° MAX SCREW ANGLULATION



CORTICAL CANCELLOUS SCREWS



ONE STEP LOCKING MECHANISM



INTEGRATED PLATE & CAGE




PRECISION SPINE®

Discover the Difference






DAKOTA ACDF™ SYSTEM



Footprints 14mm x 12mm & 16mm x 14mm
Lordosis 7°
Heights 6 - 12mm, 1mm increments



Integrated Plate and Cage Interbody



CORTICAL CANCELLOUS SCREWS







Self-Drilling/Self-Tapping
3.5mm & 4.0mm



Self-Tapping
4.0mm

Lengths

-  12mm
-  14mm
-  15mm
-  17mm

Precision Spine, Inc.

2050 Executive Drive, Pearl, MS 39208
 Customer Service: 1.888.241.4773
 Phone: 601.420.4244
 Toll Free: 877.780.4370
 Fax: 601.420.5501
www.precisionspineinc.com

