

PRECISION SPINE FACET SCREW SYSTEM

ALTERNATIVE TO PEDICLE SCREW FIXATION; STABILIZES SPINE, MINIMAL INVASION

The Facet Screw System is intended to stabilize the spine, by trans-facet fixation, as an aid to fusion. The system is composed of Facet Screws in two diameter sizes and six lengths.

The implant and associated instruments are designed to provide stabilization to the spine with minimal invasion, and offers an alternative to pedicle screw fixation.

The Facet Screw System is indicated for posterior surgical treatment with or without bone graft and can be used for single or multiple levels from L1 to S1 (inclusive) for spondylolisthesis, spondylolysis, pseudarthrosis or failed previous fusion which is symptomatic, degenerative disc disease and/or degenerative disease of the facets with instability.

System features


- Facet Screws have a self-tapping flute feature at the distal end
- The two differing threads along the length of the screw vary in aggressiveness and pitch to prevent slippage
- Screws are color coded based on diameter size for easy recognition
- Screws and instruments are designed to provide stabilization with minimal invasion
- All components of the Facet Screw system are made of medical grade titanium or titanium alloy and are supplied 'non-sterile'



FACET SCREW SYSTEM

Specifications

- 4.5mm and 5.0mm Diameters
- 25, 30, 35, 40, 45 and 50mm Lengths
- Cannulated and Non-cannulated Options



PART NO.	DESCRIPTION
35-9007	K-Wire, 1.4mm diameter x 23.5" long, trocar tip
35-9011	K-Wire, 1.4mm diameter x 23.5" long, partially threaded tip
35-9012	Introducer with handle
35-9013	Handle, Palm, Cannulated, AO connection
35-9014-8	8mm Dilator
35-9014-12	12mm Dilator
35-9015	Drill, AO connection
35-9016	Tap, AO connection
35-9017	Screw Inserter, AO connection
09-9027	K-Wire Puller

PART NUMBER	DESCRIPTION
35-S4525-K	4.5mm x 25mm Long Cannulated Screw
35-S4530-K	4.5mm x 30mm Long Cannulated Screw
35-S4535-K	4.5mm x 35mm Long Cannulated Screw
35-S4540-K	4.5mm x 40mm Long Cannulated Screw
35-S4545-K	4.5mm x 45mm Long Cannulated Screw
35-S4550-K	4.5mm x 50mm Long Cannulated Screw
35-S5025-K	5.0mm x 25mm Long Cannulated Screw
35-S5030-K	5.0mm x 30mm Long Cannulated Screw
35-S5035-K	5.0mm x 35mm Long Cannulated Screw
35-S5040-K	5.0mm x 40mm Long Cannulated Screw
35-S5045-K	5.0mm x 45mm Long Cannulated Screw
35-S5050-K	5.0mm x 50mm Long Cannulated Screw
35-S4525-X	4.5mm x 25mm Long Non-Cannulated Screw
35-S4530-X	4.5mm x 30mm Long Non-Cannulated Screw
35-S4535-X	4.5mm x 35mm Long Non-Cannulated Screw
35-S4540-X	4.5mm x 40mm Long Non-Cannulated Screw
35-S4545-X	4.5mm x 45mm Long Non-Cannulated Screw
35-S4550-X	4.5mm x 50mm Long Non-Cannulated Screw
35-S5025-X	5.0mm x 25mm Long Non-Cannulated Screw
35-S5030-X	5.0mm x 30mm Long Non-Cannulated Screw
35-S5035-X	5.0mm x 35mm Long Non-Cannulated Screw
35-S5040-X	5.0mm x 40mm Long Non-Cannulated Screw
35-S5045-X	5.0mm x 45mm Long Non-Cannulated Screw
35-S5050-X	5.0mm x 50mm Long Non-Cannulated Screw

Precision Spine, Inc.
 2050 Executive Drive, Pearl, MS 39208
 Customer Service: 1.888.241.4773
 Phone: 601.420.4244
 Toll Free: 877.780.4370
 Fax: 601.420.5501
www.precisionspineinc.com

