

# Ovation<sup>®</sup> Hip System



**STANDARD AND EXTENDED SIZES - OVATION AND OVATION TRIBUTE**

SIZE	OVATION STEM LENGTH	TRIBUTE STEM LENGTH	NECK ANGLE	HORIZONTAL OFFSET (STD)	HORIZONTAL OFFSET (EXT)	VERTICAL OFFSET	PROXIMAL WIDTH (M/L)	OVATION DISTAL WIDTH (M/L)	TRIBUTE DISTAL WIDTH (M/L)
1	98	74	132°	34	39	29	26	5	8
2	106	77	132°	34	39	29.5	27	6	9
3	108	80	132°	35	40	30	28	6	9
4	110	83	132°	37	43	32	29	7	10
5	114	86	132°	38	44	32	30	8	11
6	118	89	132°	38	44	33	31	9	12
7	122	92	132°	42	49	34	32	10	13
8	126	95	132°	43	50	34	34	11	14
8.5	—	96.5	132°	43	50	34.5	35	—	15
9	130	98	132°	44	51	34.5	36	13	16
9.5	—	99.5	132°	44	51	35	37	—	17
10	134	101	132°	47	55	36	38	16	18
10.5	—	102.5	132°	47	55	36	39	—	19
11	138	104	132°	48	56	36	40	18	20
11.5	—	105.5	132°	48	56	36	41	—	21
12	142	—	132°	49	57	36	42	20	—

Ovation and Ovation Tribute stem dimensions and offsets (all dimensions are in mm)



## Ovation® Hip System

- M/L Tapered stem design provides excellent rotational stability and promotes 3-point fixation
- Accommodates multiple surgical approaches, including Direct Anterior
- Minimal Broach-Only instrumentation simplifies surgical procedure
- Standard and reduced lengths to address proximal/distal mismatch



Ortho Development Corporation  
odev.com

12187 S Business Park Drive  
Draper, Utah 84020  
801-553-9991/fax 801-553-9993