


PRECISION SPINE  
**SLIMPLICITY**<sup>®</sup>  
ANTERIOR CERVICAL PLATING SYSTEM

**SIMPLE & VERSATILE**


The Slimplicity<sup>®</sup> Anterior Cervical Plate offers one of the slimmest plates available with an easy-to-use locking mechanism that facilitates visual locking confirmation.

Large graft windows have been incorporated to provide an unimpeded graft site and end plate visualization.

The large array of variable and fixed screw options accommodates semi-constrained, constrained and hybrid philosophies, while the plates are pre-contoured to address varying patient anatomy in multi-level constructs.



# SLIMPLICITY® ANTERIOR CERVICAL PLATING SYSTEM



**Low Profile Plate (2mm)**  
designed to minimize tissue disruption and post-op discomfort

**Pre-Contoured**  
to address patient anatomy



**Central Graft Windows**  
facilitate unimpeded graft visualization

**One Step Locking**  
facilitates secure application

## IMPLANT SIZING

### Plate Sizing

Lengths\*

- 1 Level - 7 sizes, 20-32mm (2mm increments)
- 2 Level - 8 sizes, 37-55mm (3mm increments)
- 3 Level - 8 sizes, 54-75mm (3mm increments)
- 4 level - 6 sizes, 69-89mm (4mm increments)


\*Measured end-to-end

Width 17mm

Waist 14mm


### Fixed & Variable Screw Sizes

- 4.0 and 4.5mm Diameters
- 12-18mm Lengths  
(2mm increments)



Single Tray Layout minimizes case complexity

## PRECISION SPINE SHURFIT® CERVICAL INTERBODY CAGES



- **Large Contact Area**  
designed to optimize vertebral body support and minimize risk of subsidence
- **Large Graft Area**  
enables excellent biological coverage
- **Trapezoidal Design**  
allows for proper anterior placement
- **Tantalum Markers**  
strategically placed to aid interbody positioning

### SIZING

- 14mm Wide x 12mm Deep
- 16mm Wide x 14mm Deep
- Heights - 5-12mm (1mm increments)
- 5° & 10° Lordotic Options

**Precision Spine, Inc.**  
2050 Executive Drive, Pearl, MS 39208  
Customer Service: 1.888.241.4773  
Phone: 601.420.4244  
Toll Free: 877.780.4370  
Fax: 601.420.5501  
www.precisionspineinc.com

